



Oxygenated Facelift System

... A Non-Invasive Tool to Beautify Skin

The OxyLift uses a safe, effective high frequency current to tone the skin, which not only rejuvenates your appearance but actually makes your skin healthier! This proven technology has been used for decades both in North America and Europe for therapeutic devices to help heal injuries and esthetic salon treatments to revitalize and beautify the face and neck.



The OxyLift stimulates the skin's layers with a rapid massaging action to increase blood circulation and cell growth, promote healing of blemishes, and even stimulate new collagen. Decreased collagen production is a major cause of aging in the skin so stimulating new growth prevents sagging and wrinkles and provides a youthful, vibrant glow. It reduces puffiness, decreases pore size and strengthens the skin. The OxyLift hand-held device offers a "non-surgical facelift" because it creates a temporary lifting of facial muscles, much like a workout at the gym increases tone of the body's muscles. And, like a workout, when done regularly, the effects are cumulative.

The high frequency current used in the device produces oxygen (ozone), which not only increases circulation to enliven a sallow skin tone, it kills bacteria so that skin conditions including blemishes and acne are reduced. Its therapeutic action, which increases new cell growth, is ideal for treating burns — it has been used for this purpose by specialists for decades. The increased blood circulation achieved by OxyLift will also help to gently dissipate broken capillaries and sebaceous cysts on facial skin. It has even been shown to stimulate hair growth! Cellulite and stretch marks on other areas of the body have also been reduced by OxyLift's stimulating action.

Take the specialized salon treatment home! With the handheld OxyLift device, you can create a salon experience with all of the benefits in the comfort of your bedroom.

To Use the Device ...



To use the device, introduce your facial skin to its action slowly, using it gently on the lowest setting twice a day for 2 to 3 minutes per treatment. As skin becomes accustomed to the high frequency current's action, increase the intensity. However, it is not necessary to increase the treatment time over 3 to 5 minutes per session. Use this tool every day for short intervals for maximum effect.

Begin by gently cleansing the skin and then applying your moisturizer, and then proceed with the treatment. The OxyLift will actually help your skincare products to penetrate the skin's layers more thoroughly, therefore increasing their effectiveness.

If you have oily skin, only apply moisturizer after using your OxyLift, not before. If you have normal skin and it becomes dry during the treatment, it may be due to not using enough moisturizer or from using the tool for too long within one session. Caution should be taken with rosacea skin conditions due to the broken capillaries which are close to the skin's surface. Use the lowest setting and work around these sensitive areas. Please refer to the Usage Guide and Instructional DVD that are included in your OxyLift kit for exercises and techniques that target specific areas and concerns.

Independent **Clinical Studies** Showed:

- **78%** Improvement in 5 days for acne skin conditions
- **56%** Decrease in pore size within 14 days
- **33%** Increase of collagen production in 30 days
- **43%** Decrease in puffy eyes (due to fluid retention)
- **48%** Reduction in wrinkle depth (12 weeks)
- **92%** of subjects experienced immediate lift response

Oxylift™ is a clinically proven system with corrective results for:

- **Acne**
- **Enlarged pores**
- **Puffy eyes**
- **Saggy skin**
- **Wrinkles**
- **Thinning Hair**
- **Cellulite**
- **Case Studies**
- **Statistics**



More Benefits of Using Oxylift™

- Although you receive instant results, the benefits are accumulative so you will have long term positive effects.
- The cells are rejuvenated and regenerated.
- Collagen is what produces stronger skin. Stronger skin resists wrinkles and sagging.
- The enriched oxygen doesn't allow for bacteria which is why **Oxylift™** is so effective for acne.
- O3 neutralizes the pollutants out of the skin which are causing our skin to age 10 years faster.
- O3 kills the bacteria and changes the molecular structure of our skin care products so they actually penetrate the skin better and can work from the inside out while **Oxylift™** works from the outside in.
- **Oxylift™** eliminates the doughy, tired look through a retexturing effect giving your skin the porcelain look.
- Men's skin is 1/3 thicker which is why a man of 50 will inevitably look younger than a woman of 50. Thin skin can't defend itself from bacteria, pollutants and aging. Many anti-aging products and techniques are actually thinning the skin.
- The **Oxylift™** strengthens the skin and as it has memory, it becomes stronger and therefore more resistant to wrinkles every day.

Oxylift™ is not a luxury, it's a necessity. Are you ready to try your **Oxylift™** machine today? Are you ready to look younger and feel better?